

兔抗 CLDN6(Phospho-Tyr219) 多克隆抗体

中文名称：兔抗 CLDN6(Phospho-Tyr219) 多克隆抗体

英文名称：Anti-CLDN6(Phospho-Tyr219) rabbit polyclonal antibody

别 名：CLDN6

相关类别：一抗

储 存：冷冻（-20℃）避光

宿 主：Rabbit

抗 原：CLDN6(Phospho-Tyr219)

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:

Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet, with complementary grooves in the inwardly facing extracytoplasmic leaflet. This gene encodes a component of tight junction strands, which is a member of the claudin family. The protein is an integral membrane protein and is one of the entry cofactors for hepatitis C virus.

	rus. The gene methylation may be involved in esophageal tumorigenesis. This gene is adjacent to another family member CLDN9 on chromosome 16.[provided by RefSeq, Aug 2010]
Applications:	WB
Name of antibody:	CLDN6(Phospho-Tyr219)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine 219 (T-K-N-Y(p)-V) derived from Human Claudin 6.
Full name:	claudin 6
Synonyms :	CLDN6
SwissProt:	P56747
WB Predicted band size:	23 kDa
WB Positive control:	RAW264.7 cells lysates treated with UV
WB Recommended dilution:	500-1000

