

兔抗 CLDN7 多克隆抗体

中文名称：兔抗 CLDN7 多克隆抗体

英文名称：Anti-CLDN7 rabbit polyclonal antibody

别名：CLDN-7, CEPTRL2, CPETRL2, Hs.84359, claudin-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CLDN7

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. Differential expression of this gene has been observed in different types of malignancies, including breast cancer, ovarian cancer, hepatocellular carcinomas, urinary tumors, prostate cancer, lung cancer, head and neck cancers, thymo-

	id carcinomas, etc.. Alternatively spliced transcript variants encoding different isoforms have been found.
Applications:	ELISA, WB
Name of antibody:	CLDN7
Immunogen:	Synthetic peptide of humanCLDN7
Full name:	claudin 7
Synonyms :	CLDN-7, CEPTRL2, CPETRL2, Hs.84359, claudin-1
SwissProt:	O95471
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	23 kDa
WB Positive control:	Mouse liver tissue
WB Recommended dilution:	200-1000

