

兔抗 CLDN25 多克隆抗体

中文名称: 兔抗 CLDN25 多克隆抗体

英文名称: Anti-CLDN25 rabbit polyclonal antibody

储 存: 冷冻(-20℃)

抗 原: CLDN25

宿 主: Rabbit

反应种属: Human

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the claudin family. Cla udins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epit helial or endothelial cell sheets, and also play critical rolles in maintaining cell polarity and signal transductions.
Applications:	ELISA, WB, IHC
Name of antibody:	CLDN25
Immunogen:	Synthetic peptide of human CLDN25
Full name:	claudin 25
SwissProt:	C9JDP6
ELISA Recommended dilution:	2000-5000



IHC positive control:	Human ovarian cancer and Human esophagus cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	25 kDa
WB Positive control:	Human normal lung tissue
WB Recommended dilution:	500-2000





