

兔抗 CLDND2 多克隆抗体

中文名称: 兔抗 CLDND2 多克隆抗体

英文名称: Anti-CLDND2 rabbit polyclonal antibody

抗 原: CLDND2

储 存: 冷冻 (-20℃)

宿 主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	CLDND2 (claudin domain containing 2) is a 167 amino acid multi-pass membrane protein that belongs to the PMP-22/EMP/MP20 family and is encoded by a gene that maps to human chromosome 19q13.33. Consisting of around 63 million bases with over 1,400 genes, chromosome 19 makes up over 2% of human genomic DNA. Chromosome 19 includes a diversity of interesting genes and is recognized for having the greatest gene density of the human chromosomes.
Applications:	ELISA, WB
Name of antibody:	CLDND2
Immunogen:	Synthetic peptide of human CLDND2

Full name:	claudin domain containing 2
SwissProt:	Q8NHS1
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	18 kDa
WB Positive control:	TM4, PC3, 293T and hela cell
WB Recommended dilution:	200-1000



