

兔抗 CLEC1B 多克隆抗体

中文名称：兔抗 CLEC1B 多克隆抗体

英文名称：Anti-CLEC1B rabbit polyclonal antibody

别名：C-type lectin domain family 1 member B; CLEC2; CLEC2B; PRO1384; QDED721; 1810061I13Rik

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CLEC1B

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Natural killer (NK) cells express multiple calcium-dependent (C-type) lectin-like receptors, such as CD94 (KLRC1; MIM 602894) and NKG2D (KLRC4; MIM 602893), that interact with major histocompatibility complex class I molecules and either inhibit or activate cytotoxicity and cytokine secretion. CLEC2 is a C-type lectin-like receptor expressed in myeloid cells and NK cells (Colonna et al., 2000 [PubMed 10671229]).
Applications:	ELISA, IHC
Name of antibody:	CLEC1B

Immunogen:	Fusion protein of human CLEC1B
Full name:	C-type lectin domain family 1 member B
Synonyms:	CLEC2; CLEC2B; PRO1384; QDED721; 1810061I13Rik
SwissProt:	Q9P126
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human liver cancer
IHC Recommend dilution:	50-200



