

兔抗 CMPK1 多克隆抗体

中文名称：兔抗 CMPK1 多克隆抗体

英文名称：Anti-CMPK1 rabbit polyclonal antibody

别名：cytidine/uridine monophosphate kinase 1; CK; CMK; UMK; CMPK; UMPK; UMP-CMPK

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CMPK1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes one of the enzymes required for cellular nucleic acid biosynthesis. This enzyme catalyzes the transfer of a phosphate group from ATP to CMP, UMP, or dCMP, to form the corresponding diphosphate nucleotide. Alternate splicing results in both coding and non-coding transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	CMPK1
Immunogen:	Full length fusion protein
Full name:	cytidine/uridine monophosphate kinase 1

Synonyms:	CK; CMK; UMK; CMPK; UMPK; UMP-CMPK
SwissProt:	P30085
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human esophagus cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	22 kDa
WB Positive control:	HT-29, Hela, HEPG2, A549, Jurkat and 231 cell lysate
WB Recommended dilution:	200-1000



