

## CA1 抗原（重组蛋白）

中文名称：CA1 抗原（重组蛋白）

英文名称：CA1 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	Carbonic anhydrase 1
<b>Synonyms:</b>	CAB, CA-I, Car1
<b>Swissprot:</b>	P00915
<b>Gene Accession:</b>	BC027890
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. They show extensive diversity in tissue distribution and in their subcellular localization. CA1 is closely linked to CA2 and CA3 genes on chromosome 8, and it encodes a cytosolic protein which is found at the highest level in erythrocytes. Variants of this gene have been described in some populations. Multiple alternatively spliced

ed variants, encoding the same protein, have been identified. Transcript variants of CA1 utilizing alternative polyA\_sites have been described in literature.