

兔抗 ARRDC2 多克隆抗体

中文名称：兔抗 ARRDC2 多克隆抗体

英文名称： Anti-ARRDC2 rabbit polyclonal antibody

别 名： PP2703; CLONE24945

抗 原： ARRDC2

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

相关类别： 一抗

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Belongs to the arrestin family. 2 isoforms of the human protein are produced by alternative splicing.
Applications:	IHC
Name of antibody:	ARRDC2
Immunogen:	Synthesized peptide derived from internal of human ARRDC2.
Full name:	arrestin domain containing 2
Synonyms :	PP2703; CLONE24945
SwissProt:	Q8TBH0
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100

