

## 兔抗 ARSA 多克隆抗体

中文名称：兔抗 ARSA 多克隆抗体

英文名称：Anti-ARSA rabbit polyclonal antibody

别名：MLD

相关类别：一抗

抗原：ARSA

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene hydrolyzes cerebroside sulfate to cerebroside and sulfate. Defects in this gene lead to metachromatic leucodystrophy (MLD), a progressive demyelination disease which results in a variety of neurological symptoms and ultimately death. Alternatively spliced transcript variants have been described for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ARSA
<b>Immunogen:</b>	Fusion protein of human ARSA
<b>Full name:</b>	arylsulfatase A

<b>Synonyms :</b>	MLD
<b>SwissProt:</b>	P15289
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human brain
<b>IHC Recommend dilution:</b>	50-200

