

兔抗 ARSD 多克隆抗体

中文名称：兔抗 ARSD 多克隆抗体

英文名称：Anti-ARSD rabbit polyclonal antibody

别名：ASD

抗原：ARSD

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The protein encoded by this gene is a member of the sulfatase family. Sulfatases are essential for the correct composition of bone and cartilage matrix. The encoded protein is posttranslationally glycosylated and localized to the lysosome. This gene is located within a cluster of similar arylsulfatase genes on chromosome X. A related pseudogene has been identified in the pseudoautosomal region of chromosome Y.
Applications:	IHC, IF
Name of antibody:	ARSD
Immunogen:	Synthesized peptide derived from internal of human ARSD.
Full name:	arylsulfatase D

Synonyms :	ASD
SwissProt:	P51689
IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100
IF positive control:	MCF-7 cells
IF Recommend dilution:	100-500



