

## CEP70 抗原（重组蛋白）

中文名称： CEP70 抗原（重组蛋白）

英文名称： CEP70 Antigen (Recombinant Protein)

别名： centrosomal protein 70; BITE

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-212 amino acids of human CEP70

技术规格

<b>Full name:</b>	centrosomal protein 70
<b>Synonyms:</b>	BITE
<b>Swissprot:</b>	Q8NHQ1
<b>Gene Accession:</b>	BC016050
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Plays a role in the organization of both preexisting and nascent microtubules in interphase cells. During mitosis, required for the organization and orientation of the mitotic spindle.