

兔抗 CPNE3 多克隆抗体

中文名称：兔抗 CPNE3 多克隆抗体

英文名称：Anti-CPNE3 rabbit polyclonal antibody

别名：copine 3 CPN3; PRO1071

抗原：CPNE3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene encodes a protein which contains two type II C2 domains in the amino-terminus and an A domain-like sequence in the carboxy-terminus. The A domain mediates interactions between integrins and extracellular ligands.
Applications:	ELISA, IHC
Name of antibody:	CPNE3
Immunogen:	Fusion protein of human CPNE3
Full name:	copine 3

Synonyms:	CPN3; PRO1071
SwissProt:	O75131
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	30-150



