

兔抗 CPSF1 多克隆抗体

中文名称：兔抗 CPSF1 多克隆抗体

英文名称：Anti-CPSF1 rabbit polyclonal antibody

别名：cleavage and polyadenylation specific factor 1 CPSF160; P/cl.18; HSU37012

相关类别：一抗

储存：冷冻（-20℃）

抗原：CPSF1

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Cleavage and polyadenylation specificity factor (CPSF) is a multisubunit complex that plays a central role in 3-prime processing of pre-mRNAs. CPSF recognizes the AAUAAA signal in the pre-mRNA and interacts with other proteins to facilitate both RNA cleavage and poly(A) synthesis. CPSF1 is the largest subunit of the CPSF complex (Murthy and Manley, 1995 [PubMed 7590244]).
Applications:	ELISA, IHC
Name of antibody:	CPSF1
Immunogen:	Fusion protein of human CPSF1

Full name:	cleavage and polyadenylation specific factor 1
Synonyms:	CPSF160; P/ci.18; HSU37012
SwissProt:	Q10570
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	20-100

