

兔抗 ATG12 多克隆抗体

中文名称：兔抗 ATG12 多克隆抗体

英文名称：Anti-ATG12 rabbit polyclonal antibody

别名：APG12; FBR93; APG12L; HAPG12

相关类别：一抗

抗原：ATG12

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Autophagy is a process of bulk protein degradation in which cytoplasmic components, including organelles, are enclosed in double-membrane structures called autophagosomes and delivered to lysosomes or vacuoles for degradation. ATG12 is the human homolog of a yeast protein involved in autophagy.
Applications:	ELISA, IHC
Name of antibody:	ATG12
Immunogen:	Fusion protein of human ATG12
Full name:	autophagy related 12
Synonyms :	APG12; FBR93; APG12L; HAPG12

SwissProt:	O94817
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human gastric cancer
IHC Recommend dilution:	50-200

