

DACT3 抗原（重组蛋白）

中文名称：DACT3 抗原（重组蛋白）

英文名称：DACT3 Antigen (Recombinant Protein)

别名：dishevelled binding antagonist of beta catenin 3; RRR1; DAPPER3

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1-167 amino acids of human DACT3

技术规格

Full name:	dishevelled binding antagonist of beta catenin 3
Synonyms:	RRR1; DAPPER3
Swissprot:	Q96B18
Gene Accession:	BC034052
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in regulation of intracellular signaling pathways during development. Specifically thought to play a role in canonical and/or non-canonical Wnt signaling pathways through interaction with DSH (Dishevelled) family proteins.