

## DHRS1 抗原（重组蛋白）

中文名称：DHRS1 抗原（重组蛋白）

英文名称：DHRS1 Antigen (Recombinant Protein)

别名：dehydrogenase/reductase 1; SDR19C1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 2-313 amino acids of human DHRS1

技术规格

<b>Full name:</b>	dehydrogenase/reductase 1
<b>Synonyms:</b>	SDR19C1
<b>Swissprot:</b>	Q96LJ7
<b>Gene Accession:</b>	BC015943
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded enzyme contains a conserved catalytic domain and likely functions as an oxidoreductase. Multiple alternatively spliced variants, encoding the same protein, have been identified.