

## 兔抗 CTAG1A/B 多克隆抗体

- 中文名称: 兔抗 CTAG1A/B 多克隆抗体
- 英文名称: Anti-CTAG1A/B rabbit polyclonal antibody
- 别 名: CTAG, CTAG1, ESO1, LAGE2, LAGE2A/CTAG1B
- 相关类别: 一抗
- 储 存: 冷冻 (-20℃)
- 宿 主: Rabbit
- 抗 原: CTAG1A/B
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

## 技术规格



Name of antibody:	CTAG1A/B
Immunogen:	Synthetic peptide of human CTAG1A/B
Full name:	Cancer/testis antigen 1
Synonyms:	CTAG, CTAG1, ESO1, LAGE2, LAGE2A/CTAG1B
SwissProt:	P78358/P78358
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	25-100

