

## 兔抗 BAG4 多克隆抗体

中文名称：兔抗 BAG4 多克隆抗体

英文名称：Anti-BAG4 rabbit ployclonal antibody

别名：SODD; BAG-4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BAG4

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the BAG1-related protein family. BAG1 is an anti-apoptotic protein that functions through interactions with a variety of cell apoptosis and growth related proteins including BCL-2, Raf-protein kinase, steroid hormone receptors, growth factor receptors and members of the heat shock protein 70 kDa family. This protein contains a BAG domain near the C-terminus, which could bind and inhibit the chaperone activity of Hsc70/Hsp70. This protein was found to be associated with the death domain of tumor necrosis factor receptor type 1 (T
--------------------	--

	NF-R1) and death receptor-3 (DR3), and thereby negatively regulates downstream cell death signaling.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	BAG4
<b>Immunogen:</b>	Synthetic peptide of human BAG4
<b>Full name:</b>	BCL2-associated athanogene 4
<b>Synonyms :</b>	SODD; BAG-4
<b>SwissProt:</b>	O95429
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	Human normal stomach tissue
<b>WB Recommended dilution:</b>	200-1000

