

## HSP90AA1 抗原（重组蛋白）

中文名称：HSP90AA1 抗原（重组蛋白）

英文名称：HSP90AA1 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 250 amino acids of human HSP90AA1

技术规格

<b>Full name:</b>	heat shock protein 90kDa alpha (cytosolic), class A member 1
<b>Synonyms:</b>	EL52; HSPN; LAP2; HSP86; HSPC1; HSPCA; Hsp89; Hsp90; LAP-2; HSP89A; HSP90A; HSP90N; HSPCAL1; HSPCAL4
<b>Swissprot:</b>	P07900
<b>Gene Accession:</b>	BC121062
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene