

兔抗 BBOX1 多克隆抗体

中文名称：兔抗 BBOX1 多克隆抗体

英文名称：Anti-BBOX1 rabbit polyclonal antibody

别名：gamma-butyrobetaine hydroxylase 1; BBH; BBOX; G-BBH; gamma-BBH

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BBOX1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes gamma butyrobetaine hydroxylase which catalyzes the formation of L-carnitine from gamma-butyrobetaine, the last step in the L-carnitine biosynthetic pathway. Carnitine is essential for the transport of activated fatty acids across the mitochondrial membrane during mitochondrial beta-oxidation.
Applications:	ELISA, WB, IHC
Name of antibody:	BBOX1
Immunogen:	Fusion protein of human BBOX1
Full name:	gamma-butyrobetaine hydroxylase 1
Synonyms:	BBH; BBOX; G-BBH; gamma-BBH

SwissProt:	O75936
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	45 kDa
WB Positive control:	Human kidney tissue lysate
WB Recommended dilution:	500-2000



