

兔抗 BBS9 多克隆抗体

中文名称：兔抗 BBS9 多克隆抗体

英文名称：Anti-BBS9 rabbit polyclonal antibody

别名：B1; D1; C18; PTHB1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BBS9

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is downregulated by parathyroid hormone in osteoblastic cells, and therefore, is thought to be involved in parathyroid hormone action in bones. The exact function of this gene has not yet been determined. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, IHC
Name of antibody:	BBS9
Immunogen:	Fusion protein of human BBS9
Full name:	Bardet-Biedl syndrome 9
Synonyms :	B1; D1; C18; PTHB1

SwissProt:	Q3SYG4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	20-100

