

## IL27 抗原（重组蛋白）

中文名称：IL27 抗原（重组蛋白）

英文名称： IL27 Antigen (Recombinant Protein)

别 名： p28; IL30; IL-27; IL27A; IL-27A; IL27p28

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein
----------------------------

技术规格

<b>Full name:</b>	interleukin 27
<b>Synonyms:</b>	p28; IL30; IL-27; IL27A; IL-27A; IL27p28
<b>Swissprot:</b>	Q8NEV9
<b>Gene Accession:</b>	BC062422
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is one of the subunits of a heterodimeric cytokine complex. This protein is related to interleukin 12A (IL12A). It interacts with Epstein-Barr virus induced gene 3 (EBI3), a protein similar to interleukin 12B (IL12B), and forms a complex that has been shown to drive rapid expansion of naive but not memory CD4(+) T cells. The complex is also found to synergize strongly with interleukin 12 to trigger interferon gamma (IFNG) production of naive CD4(+) T cells.

	The biological effects of this cytokine are mediated by the class I cytokine receptor (WSX1/TCRR)
--	---