

兔抗 BCO2 多克隆抗体

中文名称：兔抗 BCO2 多克隆抗体

英文名称：Anti-BCO2 rabbit polyclonal antibody

别名：beta-carotene oxygenase 2; BCDO2; B-DIOX-II

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BCO2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an enzyme which oxidizes carotenoids such as beta-carotene during the biosynthesis of vitamin A. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	BCO2
Immunogen:	Fusion protein of human BCO2
Full name:	beta-carotene oxygenase 2
Synonyms:	BCDO2; B-DIOX-II
SwissProt:	Q9BYV7
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human brain
IHC Recommend dilution:	25-100

