

NTMT1 抗原（重组蛋白）

中文名称：NTMT1 抗原（重组蛋白）

英文名称：NTMT1 Antigen (Recombinant Protein)

别 名：NRMT; NTM1A; AD-003; HOMT1A; C9orf32; METTL11A

储 存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	N-terminal Xaa-Pro-Lys N-methyltransferase 1
Synonyms:	NRMT; NTM1A; AD-003; HOMT1A; C9orf32; METTL11A
Swissprot:	Q9BV86
Gene Accession:	BC001396
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The METTL11A gene encodes an N-terminal methyltransferase for the RAN (MIM 601179) guanine nucleotide exchange factor regulator of chromosome condensation 1 (RCC1; MIM 179710). METTL11A enzyme alpha-N-methylates other protein targets such as SET.