

## 兔抗 BRAF(Phospho-Thr753) 多克隆抗体

中文名称：兔抗 BRAF(Phospho-Thr753) 多克隆抗体

英文名称：Anti-BRAF(Phospho-Thr753) rabbit polyclonal antibody

别名：NS7; BRAF1; RAFB1; B-RAF1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：BRAF(Phospho-Thr753)

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

#### Background:

This gene encodes a protein belonging to the raf/mil family of serine/threonine protein kinases. This protein plays a role in regulating the MAP kinase/ERKs signaling pathway, which affects cell division, differentiation, and secretion. Mutations in this gene are associated with cardiofaciocutaneous syndrome, a disease characterized by heart defects, mental retardation and a distinctive facial appearance. Mutations in this gene have also been associated with various cancers, including non-Hodgkin lymphoma, colorectal cancer, malignant melanoma, thyroid carcinoma, non-small

	cell lung carcinoma, and adenocarcinoma of lung. A pseudogene, which is located on chromosome X, has been identified for this gene. [provided by RefSeq, Jul 2008]
<b>Applications:</b>	WB
<b>Name of antibody:</b>	BRAF(Phospho-Thr753)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of threonine 753 (P-K-T(p)-P-L) derived from Human B-Raf .
<b>Full name:</b>	B-Raf proto-oncogene, serine/threonine kinase
<b>Synonyms :</b>	NS7; BRAF1; RAFB1; B-RAF1
<b>SwissProt:</b>	P15056
<b>WB Predicted band size:</b>	84 kDa
<b>WB Positive control:</b>	K562 cells lysates treated with EGF
<b>WB Recommended dilution:</b>	500-1000

