

RAP1B 抗原（重组蛋白）

中文名称：RAP1B 抗原（重组蛋白）

英文名称：RAP1B Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	Ras-related protein Rap-1b
Synonyms:	K-REV, RAL1B
Swissprot:	P61224
Gene Accession:	BC000176
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the RAS-like small GTP-binding prote in superfamily. Members of this family regulate multiple cellular proce sses including cell adhesion and growth and differentiation. This prote in localizes to cellular membranes and has been shown to regulate in tegrin-mediated cell signaling. This protein also plays a role in regulat ing outside-in signaling in platelets. Alternate splicing results in multi ple transcript variants. Pseudogenes of this gene are found on chrom osomes 3, 5, 6 and 9.