

兔抗 DACT3 多克隆抗体

中文名称：兔抗 DACT3 多克隆抗体

英文名称：Anti-DACT3 rabbit polyclonal antibody

别名：dishevelled binding antagonist of beta catenin 3; RRR1; DAPPER3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DACT3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in regulation of intracellular signaling pathways during development. Specifically thought to play a role in canonical and/or non-canonical Wnt signaling pathways through interaction with Dishevelled (Dsh) family proteins.
Applications:	ELISA, IHC
Name of antibody:	DACT3
Immunogen:	Fusion protein of human DACT3
Full name:	dishevelled binding antagonist of beta catenin 3
Synonyms:	RRR1; DAPPER3
SwissProt:	Q96B18

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

