

兔抗 BRD3 多克隆抗体

中文名称：兔抗 BRD3 多克隆抗体

英文名称：Anti-BRD3 rabbit polyclonal antibody

别名：ORFX; RING3L

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BRD3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene was identified based on its homology to the gene encoding the RING3 protein, a serine/threonine kinase. The gene localizes to 9q34, a region which contains several major histocompatibility complex (MHC) genes. The function of the encoded protein is not known.
Applications:	ELISA, IHC
Name of antibody:	BRD3
Immunogen:	Synthetic peptide of human BRD3
Full name:	bromodomain containing 3
Synonyms:	ORFX; RING3L
SwissProt:	Q15059

ELISA Recommended dilution :	5000-10000
IHC positive control:	Human tonsil and human colorectal cancer
IHC Recommend dilution:	10-50



