

兔抗 DAOA 多克隆抗体

中文名称：兔抗 DAOA 多克隆抗体

英文名称：Anti-DAOA rabbit polyclonal antibody

别名：LG72; SG72

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DAOA

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a protein that may function as an activator of D-amino acid oxidase, which degrades the gliotransmitter D-serine, a potent activator of N-methyl-D-aspartate (NMDA) type glutamate receptors. Studies also suggest that one encoded isoform may play a role in mitochondrial function and dendritic arborization. Polymorphisms in this gene have been implicated in susceptibility to schizophrenia and bipolar affective disorder. Alternatively spliced transcript variants encoding different isoforms have been identified.

Applications:	ELISA, IHC
Name of antibody:	DAOA
Immunogen:	Full length fusion protein(153FULL)
Full name:	D-amino acid oxidase activator
Synonyms:	LG72; SG72
SwissProt:	P59103
ELISA Recommended dilution:	5000-10000
IHC Positive control:	Human liver cancer;Human esophagus cancer
IHC Recommended dilution:	25-100



