

兔抗 DAPK1 多克隆抗体

中文名称：兔抗 DAPK1 多克隆抗体

英文名称：Anti-DAPK1 rabbit polyclonal antibody

别名：DAPK

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DAPK1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Death-associated protein kinase 1 is a positive mediator of gamma-interferon induced programmed cell death. DAPK1 encodes a structurally unique 160-kD calmodulin dependent serine-threonine kinase that carries 8 ankyrin repeats and 2 putative P-loop consensus sites. It is a tumor suppressor candidate. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	DAPK1
Immunogen:	Fusion protein of human DAPK1
Full name:	death-associated protein kinase 1

Synonyms :	DAPK
SwissProt:	P53355
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human esophagus cancer and human thyroid cancer
IHC Recommend dilution:	25-100



