

兔抗 BRINP1 多克隆抗体

中文名称：兔抗 BRINP1 多克隆抗体

英文名称：Anti-BRINP1 rabbit polyclonal antibody

别名：DBC1; FAM5A; DBCCR1

相关类别：一抗

储存：冷冻（-20℃）

抗原：BRINP1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is located within a chromosomal region that shows loss of heterozygosity in some bladder cancers. It contains a 5' CpG island that may be a frequent target of hypermethylation, and it may undergo hypermethylation-based silencing in some bladder cancers.
Applications:	ELISA, IHC
Name of antibody:	BRINP1
Immunogen:	Fusion protein of human BRINP1
Full name:	bone morphogenetic protein/retinoic acid inducible neural-specific 1
Synonyms:	DBC1; FAM5A; DBCCR1

SwissProt:	O60477
IHC positive control:	Human gastric cancer and human colon cancer
IHC Recommend dilution:	50-200

