

兔抗 BRMS1L 多克隆抗体

中文名称：兔抗 BRMS1L 多克隆抗体

英文名称：Anti-BRMS1L rabbit polyclonal antibody

别名：BRMS1

抗原：BRMS1L

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene shows sequence similarity to the human breast carcinoma metastasis suppressor (BRMS1) protein and the mammalian Sds3 (suppressor of defective silencing 3) proteins. This protein is a component of the mSin3a family of histone deacetylase complexes (HDAC).
Applications:	ELISA, IHC
Name of antibody:	BRMS1L
Immunogen:	Fusion protein of human BRMS1L
Full name:	breast cancer metastasis-suppressor 1-like
Synonyms:	BRMS1

SwissProt:	Q5PSV4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human prostate cancer
IHC Recommend dilution:	20-100



