

兔抗 BRMS1 多克隆抗体

中文名称: 兔抗 BRMS1 多克隆抗体

英文名称: Anti-BRMS1 rabbit polyclonal antibody

抗 原: BRMS1

储 存: 冷冻 (-20℃)

宿 主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene reduces the metastatic potential, but not the tumorigenicity, of human breast cancer and melanoma cell lines. The protein encoded by this gene localizes primarily to the nucleus and is a component of the mSin3a family of histone deacetylase complexes (HDAC). The protein contains two coiled-coil motifs and several imperfect leucine zipper motifs. Alternative splicing results in two transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	BRMS1
Immunogen:	Synthetic peptide of human BRMS1
Full name:	breast cancer metastasis suppressor 1

SwissProt:	Q9HCU9
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100

