

兔抗 DAPK3 多克隆抗体

中文名称：兔抗 DAPK3 多克隆抗体

英文名称：Anti-DAPK3 rabbit polyclonal antibody

别名：DLK; ZIP; ZIPK

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DAPK3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis.
Applications:	ELISA, IHC
Name of antibody:	DAPK3
Immunogen:	Synthetic peptide of human DAPK3
Full name:	death-associated protein kinase 3
Synonyms :	DLK; ZIP; ZIPK
SwissProt:	O43293
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human colorectal cancer and human tonsil
IHC Recommend dilution:	25-100



