

兔抗 BSG 多克隆抗体

中文名称: 兔抗 BSG 多克隆抗体

英文名称: Anti-BSG rabbit polyclonal antibody

别 名: M6; OK; 5F7; TCSF; CD147; EMMPRIN

相关类别: 一抗

储 存: 冷冻(-20℃)

抗原: BSG

宿 主: Rabbit

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a plasma membrane protein that is important in spermatogenesis, embryo imp lantation, neural network formation, and tumor progressio n. The encoded protein is also a member of the immuno globulin superfamily. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	BSG
Immunogen:	Fusion protein of human BSG
Full name:	basigin
Synonyms:	M6; OK; 5F7; TCSF; CD147; EMMPRIN



SwissProt:	P35613
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer and Human brain
IHC Recommend dilution:	50-200
WB Predicted band size:	42 kDa
WB Positive control:	A549 and 231 cells
WB Recommended dilution:	500-2000





