

兔抗 BTK(Ab-222) 多克隆抗体

中文名称：兔抗 BTK(Ab-222) 多克隆抗体

英文名称：Anti-BTK(Ab-222) rabbit polyclonal antibody

别名：AT; ATK; BPK; XLA; IMD1; AGMX1; PSCTK1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：BTK(Ab-222)

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The protein encoded by this gene plays a crucial role in B-cell development. Mutations in this gene cause X-linked agammaglobulinemia type 1, which is an immunodeficiency characterized by the failure to produce mature B lymphocytes, and associated with a failure of Ig heavy chain rearrangement. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	WB, IHC
Name of antibody:	BTK(Ab-222)
Immunogen:	Synthesized non-phosphopeptide derived from human

	BTK around the phosphorylation site of tyrosine 222 (A-L-Y(p)-D-Y).
Full name:	Bruton agammaglobulinemia tyrosine kinase
Synonyms :	AT; ATK; BPK; XLA; IMD1; AGMX1; PSCTK1
SwissProt:	Q06187
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	76 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution:	500-3000



