

兔抗 DCK 多克隆抗体

中文名称: 兔抗 DCK 多克隆抗体

英文名称: Anti-DCK rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻(-20℃)

抗原: DCK

宿 主: Rabbit

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Deoxycytidine kinase (DCK) is required for the phospho rylation of several deoxyribonucleosides and their nucle oside analogs. Deficiency of DCK is associated with resi stance to antiviral and anticancer chemotherapeutic age nts. Conversely, increased deoxycytidine kinase activity is associated with increased activation of these compounds to cytotoxic nucleoside triphosphate derivatives. DC K is clinically important because of its relationship to durug resistance and sensitivity.
Applications:	ELISA, WB, IHC
Name of antibody:	DCK
Immunogen:	Fusion protein of human DCK
Full name:	deoxycytidine kinase



SwissProt:	P27707
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100
WB Predicted band size:	31 kDa
WB Positive control:	Mouse heart tissue
WB Recommended dilution:	200-1000





