

S100PBP 抗原（重组蛋白）

中文名称：S100PBP 抗原（重组蛋白）

英文名称：S100PBP Antigen (Recombinant Protein)

别 名：S100PBPR

储 存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	S100P binding protein
Synonyms:	S100PBPR
Swissprot:	Q96BU1
Gene Accession:	BC015175
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein that was originally identified by its interaction with S100 calcium-binding protein P. Expression of this protein has been reported to be associated with pancreatic ductal adenocarcinoma. Alternatively spliced transcript variants have been described.