

兔抗 BZW2 多克隆抗体

中文名称：兔抗 BZW2 多克隆抗体

英文名称：Anti-BZW2 rabbit polyclonal antibody

别名：basic leucine zipper and W2 domains 2; MST017; HSPC028; MSTP017

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BZW2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in neuronal differentiation.
Applications:	ELISA, IHC
Name of antibody:	BZW2
Immunogen:	Fusion protein of human BZW2
Full name:	basic leucine zipper and W2 domains 2
Synonyms:	MST017; HSPC028; MSTP017
SwissProt:	Q9Y6E2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100

