

## SAA4 抗原（重组蛋白）

中文名称：SAA4 抗原（重组蛋白）

英文名称：SAA4 Antigen (Recombinant Protein)

别名：CSAA; C-SAA

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein
----------------------------

技术规格

<b>Full name:</b>	serum amyloid A4, constitutive
<b>Synonyms:</b>	CSAA; C-SAA
<b>Swissprot:</b>	P35542
<b>Gene Accession:</b>	BC007026
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	SAA4 is a member of the SAA family.SAA proteins are family of apolipoproteins of high density lipoprotein (HDL). They can be separated into two distinct groups. First group (SAA1, SAA2, and SAA3) consist s of acute phase reactant whose expression level increase in the blo od in a response to trauma, infection, inflammation, and neoplasia. T hese acute phase SAAs associates with HDL during inflammation and remodel the HDL particle by displacing Apo-A1. The second distinct

group consists of SAA4 and SAA5 which exist as the minor apolipoproteins on HDL, but this group of SAA constitutes more than 90% of all the SAA during homeostasis, and it is thought to play a role in the normal functioning of the HDL particle. SAA4 is a constitutively expressed protein which expressed only in humans and mice. It is connected almost completely with lipoproteins of the high density range. The physiological function of SAA4 is unknown, and its serum concentration has no association with those of other major apolipoproteins.