

兔抗 C1orf112 多克隆抗体

中文名称：兔抗 C1orf112 多克隆抗体

英文名称：Anti-C1orf112 rabbit polyclonal antibody

别名：chromosome 1 open reading frame 112

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：C1orf112

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	C1orf112
Immunogen:	Fusion protein of human C1orf112
Full name:	chromosome 1 open reading frame 112
SwissProt:	Q9NSG2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	100-300

