

兔抗 DCT 多克隆抗体

中文名称：兔抗 DCT 多克隆抗体

英文名称：Anti-DCT rabbit polyclonal antibody

别名：TRP-2; TYRP2

储存：冷冻（-20℃）

抗原：DCT

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This protein can involved in regulating eumelanin and phaeomelanin levels.The subcellular location of this protein is melanosome membrane, Single-pass type I membrane protein. It may form a multienzyme complex together with Tyrosinase and TYRP1.
Applications:	ELISA, IHC
Name of antibody:	DCT
Immunogen:	Fusion protein of human DCT
Full name:	dopachrome tautomerase
Synonyms :	TRP-2; TYRP2
SwissProt:	P40126

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and Human colon cancer
IHC Recommend dilution:	25-100

