

兔抗 DDA1 多克隆抗体

中文名称：兔抗 DDA1 多克隆抗体

英文名称：Anti-DDA1 rabbit polyclonal antibody

别名：DET1 and DDB1 associated 1 PCIA1; C19orf58

抗原：DDA1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in ubiquitination and subsequent proteasomal degradation of target proteins. Component of the DDD-E2 complexes which may provide a platform for interaction with CUL4A and WD repeat proteins.
Applications:	ELISA, IHC
Name of antibody:	DDA1
Immunogen:	Full length fusion protein
Full name:	DET1 and DDB1 associated 1
Synonyms:	PCIA1; C19orf58
SwissProt:	Q9BW61
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human brain
IHC Recommend dilution:	20-100

