

## SRGN 抗原（重组蛋白）

中文名称：SRGN 抗原（重组蛋白）

英文名称：SRGN Antigen (Recombinant Protein)

别名：PPG; PRG; PRG1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein
----------------------------

技术规格

<b>Full name:</b>	serglycin
<b>Synonyms:</b>	PPG; PRG; PRG1
<b>Swissprot:</b>	P10124
<b>Gene Accession:</b>	BC015516
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein best known as a hematopoietic cell granule proteoglycan. Proteoglycans stored in the secretory granules of many hematopoietic cells also contain a protease-resistant peptide core, which may be important for neutralizing hydrolytic enzymes. This encoded protein was found to be associated with the macromolecular complex of granzymes and perforin, which may serve as a mediator of granule-mediated apoptosis.