

SULT1B1 抗原（重组蛋白）

中文名称：SULT1B1 抗原（重组蛋白）

英文名称：SULT1B1 Antigen (Recombinant Protein)

别名：ST1B1; ST1B2; SULT1B2

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	sulfotransferase family, cytosolic, 1B, member 1
Synonyms:	ST1B1; ST1B2; SULT1B2
Swissprot:	O43704
Gene Accession:	BC010895
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. The gene structure (number and length of exons) is similar among family members. However, the total genomic length of this gene is greater than that of other SULT1 genes.