

TCEAL9 抗原（重组蛋白）

中文名称：TCEAL9 抗原（重组蛋白）

英文名称：TCEAL9 Antigen (Recombinant Protein)

别名：WBP5; WEX6

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	transcription elongation factor A like 9
Synonyms:	WBP5; WEX6
Swissprot:	Q9UHQ7
Gene Accession:	BC023544
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The globular WW domain is composed of 38 to 40 semiconserved amino acids shared by proteins of diverse functions including structural, regulatory, and signaling proteins. The domain is involved in mediating protein-protein interactions through the binding of polyproline ligands. This gene encodes a WW domain binding protein. This gene also encodes a domain with similarity to the transcription elongation factor A, SII-related family. Alternative splicing results in multiple transcr

ipt variants encoding a single isoform.