

## 兔抗 DDX11 多克隆抗体

中文名称：兔抗 DDX11 多克隆抗体

英文名称：Anti-DDX11 rabbit polyclonal antibody

别名：CHL1; KRG2; WABS; CHLR1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DDX11

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is an enzyme that possesses

	both ATPase and DNA helicase activities. This gene is a homolog of the yeast CHL1 gene, and may function to maintain chromosome transmission fidelity and genome stability. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DDX11
<b>Immunogen:</b>	Fusion protein of human DDX11
<b>Full name:</b>	DEAD/H (Asp-Glu-Ala-Asp/His) box helicase 11
<b>Synonyms :</b>	CHL1; KRG2; WABS; CHLR1
<b>SwissProt:</b>	Q96FC9
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human ovarian cancer and Human breast cancer
<b>IHC Recommend dilution:</b>	50-200

