

兔抗 CALML5 多克隆抗体

中文名称：兔抗 CALML5 多克隆抗体

英文名称：Anti-CALML5 rabbit polyclonal antibody

别名：calmodulin like 5; CLSP

相关类别：一抗

抗原：CALML5

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Background: | This gene encodes a novel calcium binding protein expressed in the epidermis and related to the calmodulin family of calcium binding proteins. Functional studies with recombinant protein demonstrate it does bind calcium and undergoes a conformational change when it does so. Abundant expression is detected only in reconstructed epidermis and is restricted to differentiating keratinocytes. In addition, it can associate with transglutaminase 3, shown to be a key enzyme in the terminal differentiation of keratinocytes. |
| Applications: | ELISA, IHC |

| | |
|------------------------------------|--------------------------------------------------|
| Name of antibody: | CALML5 |
| Immunogen: | Full length fusion protein |
| Full name: | calmodulin like 5 |
| Synonyms: | CLSP |
| SwissProt: | Q9NZT1 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human colorectal cancer and human ovarian cancer |
| IHC Recommend dilution: | 50-300 |



