

兔抗 CAMK4 多克隆抗体

中文名称：兔抗 CAMK4 多克隆抗体

英文名称：Anti-CAMK4 rabbit polyclonal antibody

别名：IV; caMK; CaMK IV; CaMK-GR

相关类别：一抗

储存：冷冻（-20℃）

抗原：CAMK4

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The product of this gene belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This enzyme is a multifunctional serine/threonine protein kinase with limited tissue distribution, that has been implicated in transcriptional regulation in lymphocytes, neurons and male germ cells.
Applications:	ELISA, IHC
Name of antibody:	CAMK4
Immunogen:	Fusion protein of human CAMK4
Full name:	calcium/calmodulin-dependent protein kinase IV
Synonyms:	IV; caMK; CaMK IV; CaMK-GR

SwissProt:	Q16566
IHC positive control:	Human gastric cancer and human breast cancer
IHC Recommend dilution:	50-200

