

兔抗 DGAT1 多克隆抗体

中文名称：兔抗 DGAT1 多克隆抗体

英文名称：Anti-DGAT1 rabbit polyclonal antibody

别名：ARAT; DGAT; ARGP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DGAT1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an multipass transmembrane protein that functions as a key metabolic enzyme. The encoded protein catalyzes the conversion of diacylglycerol and fatty acyl CoA to triacylglycerol. This enzyme can also transfer acyl CoA to retinol. Activity of this protein may be associated with obesity and other metabolic diseases.
Applications:	ELISA, IHC
Name of antibody:	DGAT1
Immunogen:	Synthetic peptide of human DGAT1
Full name:	diacylglycerol O-acyltransferase 1
Synonyms:	ARAT; DGAT; ARGP1

SwissProt:	O75907
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	25-100

